



# Arisure<sup>®</sup> Dry Spike

The Arisure<sup>®</sup> Dry Spike gives you an easy, safe way to transfer fluid into a IV bag or container while maintaining a closed system. When used with the Arisure<sup>®</sup> Closed Male Luer (CML), the device provides a secure, intuitive connection and a dry disconnection that reduces exposure to harmful drug residue. A double seal design allows for secure and safe connections with a variety of IV set spikes. The device's high flow rate facilitates the efficient transfer of medication from a bag, bottle, or other semi-rigid containers to the patient

The Arisure<sup>®</sup> Dry Spike forms part of the Arisure<sup>®</sup> Closed System Transfer System, designed for healthcare professionals to ensure the safe compounding, transport, and delivery of hazardous drugs.

## Features and Benefits

- ✓ Simple and easy to attach securely to a bag, bottle, or other semi-rigid container
- ✓ Dry disconnect when used with Arisure<sup>®</sup> CML results in reduced exposure to harmful drug residue
- ✓ Double seal design allows for secure, safe connection with a variety of spikes
- ✓ High flow rate means an easy transfer of drug or dose into the container



## Ordering Information

### IV0050

Arisure<sup>®</sup> Dry Spike  
sterile • 50 pcs/pck



Enquire  
today for  
**TRIAL**  
samples



## Technical Data

	IV0050
Material	Thermoplastic Elastomer; Acrylonitrile-Butadiene-Styrene; Low-Density Polymer; Polycarbonate; Solicone
Luer Access	Compatible with male luers; Arisure® Closed Mail Luer required for dry disconnect and closed system transfer
Properties	Dual lumen spike for drug admixture and patient delivery; Not constructed with animal content, Phthalates, PVC or DEHP; Latex-Free
Flow rate (Needle-free Valve)	>32ml/min
Flow rate (Spike Lumen)	>100ml/min
Residual volume	1.5ml
Drug Compatibility	Performance evaluated with a list of common hazardous medications*
Packaging	Individually blister packed; 50 x blisters per pack
Sterilisation	Gamma sterilised
Compliance	CE marked; 510(k) cleared by United States Food and Drug Administration
Shelf Life	3 years from date of manufacture

\*Data available on request

## The Complete Arisure® Closed Transfer System



IV0051 / IV0053  
Arisure Closed Vial Adapters



IV0058  
Arisure Closed Male Luer



IV0050  
Arisure Dry Spike

### Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

**Web:** helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

