



Scan Me for Web Page



Blunt Fill Aspirator Needles

Ideal for the precise withdrawal and introduction of fluids into deep well containers and vials, Helapet Blunt Fill Aspirator needles use a 16 gauge blunt cannula to prevent unnecessary needle-stick injury. Available in two convenient cannula lengths, each needle incorporates a needle guard for added protection during transit and safe transfer into controlled environments.

Features and Benefits

- ✓ Ideal for use with ampoules, vials and bottles
- ✓ Reduces risk of needlestick injury
- ✓ Easy penetration of vial stoppers
- ✓ Two convenient cannula length options
- ✓ Sterile
- ✓ CE marked



Ordering Information

IV3234

3.5"
sterile • 100 pcs/pck

IV3231

5"
sterile • 100 pcs/pck

Technical Data

Material	Stainless steel 16g blunt needle; Polypropylene Hub
Needle Guard	Polyvinyl Chloride
Cannula (IV3234)	3.5" Length
Cannula (IV3231)	5" Length
Properties	Latex-free
Packaging	Individually Blister Packed
Sterilisation	Ethylene Oxide (ETO)
Compliance	CE marked
Shelf Life	3 years from date of manufacture

Enquire today for
TRIAL
samples