



Enteral Feeding Consumables for safe access and patient delivery

A photograph of a two-story brick building with large windows. A Helapet sign is mounted on the wall, and a 'RECEPTION VISITOR PARKING' sign is visible. A 'WARNING' sign is also present.

Welcome to Helapet

Helapet has been established for over 40 years, supplying quality medical consumables to the NHS and private healthcare. Our products are design for provide top quality patient care, supporting practitioners in the preparation, handling and administration of medications and other substances.

We are proud to offer a dedicated selection of enteral feeding devices, supporting healthcare professional in the delivery of nutrition to both adult and paediatric patients. Our products meet the safety and non-interconnectivity criteria of ISO 80369-1 and ISO 80369-3.

OUR CUSTOMER COMMITMENT

At Helapet we pride our success and growth on our ability to put the customer first. By placing our customer at the centre of our brand, we continually seek to understand their wants and needs, enabling us to provide the best customer experience. We have adopted a highly customer-centric culture where all Helapet staff are trained and empowered to improve how we do things for the ultimate benefit of our clients.



OUR KNOW-HOW AND EXPERTISE

Helapet uses our wealth of expertise to deliver a high-quality customer experience, fostering positive, long-term relationships with our clients. Our service places accessibility, close consultancy, education and technical support at the forefront of what we do, to effectively understand and help our customers through their journey.

Email sales@helapet.co.uk to arrange an on-site visit from your dedicated territory Account Manager.




Learn more about us by watching our video

“ Always found Helapet to be extremely efficient and helpful”

District General Hospital

Enteral Syringes

1ml to 50ml syringes for enteral medication, rinsing, removal of gastric residues and for administration of milk or nutrition by gravity, bolus or syringe pumps.



Nutricair™ Enteral syringes

Nutricair™ Enteral syringes are transparent, with double graduation for 5ml, 10ml, 20ml, 50ml capacities, ENFit female tip, and ENFit female low dose tip for 1ml and 2.5ml syringes.

Product Code	Volume (ml)	Tip position	Sales Unit
NCE01-V2SE	1	Low Dose Tip, centric	100
NCE02SE	2.5		
NCE05SE	5	centric	50
NCE10SE	10		
NCE20SE	20		
NCE50SE	50		

Accessories

Caps for Nutricair™ Enteral syringes

Product Code	Description	Sales Unit
NCE221C	Snap-on cap for ENFit female tip x 1 unit	200
NCE228C	Snap-on cap for ENFit female tip x 8 units	480

MD class: IIa, CE0459. Please read the instruction for use



NCE221C

Nutricair™ Enteral cap for feeding tubes

Product Code	NHSSC NPC Code†	Description	Sales Unit
NCE301C	FAS85068	ENFit cap for feeding tube x 1 unit	200

MD class: IIa, CE0459. Please read the instruction for use



NCE301C

Nutricair™ Enteral bottle caps for drug sampling

Product Code	NHSSC NPC Code†	Description	Sales Unit
NCE20C	-	Nutricair™ bottle cap for enteral nutrition with ENFit for vial Ø 20 mm	40
NCE100C	-	Nutricair™ bottle cap for enteral nutrition with ENFit for vial Ø 10-24.5 mm	100
NCE302C	FAS85071	Nutricair™ screw cap for enteral nutrition with ENFit for vial Ø 28 mm	50

MD class: Is, CE0459. Please read the instruction for use



NCE20C

NCE100C

Nutricair™ Enteral transition adapters for enteral feeding

Product Code	NHSSC NPC Code†	Description	Sales Unit
NCE100A	-	Soft funnel with cap - ENFit female	150
NCE105A	FAS85069	ENFit male - stepped connector	100

MD class: Is, CE0459. Please read the instruction for use



NCE100A

NCE105A

Nutricair™ Enteral stopcock

Product Code	NHSSC NPC Code†	Description	Sales Unit
NCE300R	FAS85057	3-way stopcock: 2 ENFit males, 1 ENFit female	50



NCE301C

Nutricair™ Enteral Y port adapters

Product Code	NHSSC NPC Code†	Description	Sales Unit
NCE100Y	FAS85055	2 ENFit male and 1 ENFit female	100
NCE102Y	-	1 ENFit male, 1 funnel and 1 stepped connector	50

MD class: IIa, CE0459. Please read the instruction for use



NCE100Y

Nutricair™ Enteral adapter for enteral reservoir

Product Code	Description	Sales Unit
NCT514A	Cap for water bottle with ENPlus cross port	100

† NPC - National Product Code



Gravity sets



NCE9011YHN

Gravity tubing with ENPlus cross-over trocar for ENPlus flexible enteral feeding reservoir with ENFit female distal tip

Product Code	Description	Sales Unit
NCE9011N	with ENPlus + ENFit male medication port + ENFit female distal port	60
NCE9011YHN	Y access with ENPlus and water bottle adapter + ENFit male medication port + ENFit female distal port	40
NCE9012H	for water bottle	25
NCE9012N	with ENPlus + ENFit female distal port	60

MD class: Is, CE0459. Please read the instruction for use

Pump sets

ENPlus cross-table pump tubing for ENPlus flexible enteral reservoir with ENFit male medication/rinse port and ENFit female distal tip

Product Code	Pump	Description	Sales Unit
NCE4002N	Compat, NESTLE/NOVARTIS	ENPlus Y-tubing with ENFit medication/rinse port	60
NCE6002N	Vial mini 1000/2000, FRESenius VIAL		40
NCE6002YHN			60
NCE7002N	P7000, NESTLE		40
NCE7002YHN			60
NCE7802N	Flocare 800, NUTRICIA		60

MD class: IIa, CE0459. Please read the instruction for use

Nutricair™ Enteral feeding tubes - Neonatology / Paediatrics

Nutricair™ Enteral feeding tubes are available in Polyvinyl Chloride (PVC) or Polyurethane (PU) versions. They are transparent with radio-opaque border, centimetric numerical marking, atraumatic closed distal tip, 2 lateral distal openings, ENFit male connector with screwed cap and coloured tether per French size.



• Nutricair™ Enteral PVC naso-gastric feeding tubes

Product Code	NHSSC NPC Code†	Ch/Fr	Cap tether colour-code	Length (cm)	Sales Unit
NCE04P	FAS85065	4	Purple	40	40
NCE104P	-			125	
NCE05P	FAS85053	5	Light grey	50	
NCE105P	-			125	
NCE06P	FAS85075	6	Light green	50	
NCE106P	-			125	
NCE08P	FAS85059	8	Blue	50	
NCE108P	-			125	
NCE10P	-	10	Dark grey	50	
NCE110P	FAS85064			125	

• Nutricair™ Enteral PU naso-gastric feeding tubes

Product Code	NHSSC NPC Code†	Ch/Fr	Cap tether colour-code	Length (cm)	Sales Unit
NCE04PU	FAS85058	4	Purple	40	40
NCE104PU	-			125	
NCE05PU	FAS85067	5	Light grey	50	
NCE105PU	-			125	
NCE06PU	-	6	Light green	50	
NCE106PU	-			125	
NCE306PU	-	6	Light green	110	
NCE08PU	-			50	
NCE108PU	-			125	
NCE308PU	-			110	
NCE10PU	-	10	Dark grey	50	
NCE110PU	-			125	
NCE310PU	-			110	
NCE312PU	-			110	

MD class: IIa, CE0459

† NPC - National Product Code

Nutricair™ Enteral feeding tubes - Adults

Nutricair™ Enteral feeding tubes in Polyurethane (PU) are transparent with radio-opaque line, centimetre marking and numerical marking every 5cm (from 20 to 100 cm), atraumatic closed distal tip, 4 lateral distal openings, ENFit male connector with screwed cap and coloured tether.



NCE506AU

• Nutricair™ Enteral naso-gastric feeding tubes

Product Code	NHSSC NPC Code†	Ch/Fr	Cap tether color-code	Length (cm)	Sales Unit
NCE106AU	FAS85062	6	Light green	110	20
NCE506AU	-			130	
NCE108AU	FAS85073	8	Blue	110	
NCE508AU	-			130	
NCE110AU	FAS85063	10	Dark grey	110	
NCE510AU	-			130	
NCE112AU	FAS85052	12	White	110	
NCE512AU	-			130	
NCE114AU	FAS85070	14	Dark green	110	
NCE514AU	-			130	
NCE116AU	FAS85056	16	Orange	110	
NCE516AU	-			130	
NCE118AU	-	18	Red	110	
NCE518AU	-			130	
NCE120AU	-	20	Yellow	110	
NCE520AU	-			130	

• Nutricair™ Enteral naso-gastric feeding tubes with guide-wire

Product Code	NHSSC NPC Code†	Ch/Fr	Cap tether colour-code	Length (cm)	Sales Unit
NCE106G	FAS85054	6	Light green	40	40
NCE506G	-			125	
NCE108G	FAS85060	8	Blue	50	
NCE508G	-			125	
NCE110G	FAS85066	10	Dark grey	50	
NCE510G	-			125	
NCE112G	FAS85061	12	White	50	
NCE512G	-			125	
NCE114G	FAS85074	14	Dark green	50	
NCE514G	-			125	
NCE116G	FAS85072	16	Orange	50	
NCE516G	-			125	

† NPC - National Product Code

Gastric adult suction tubes

The Salem-type gastric suction tubes, double channel, made of Polyurethane (PU) are transparent with radio-opaque line, centimetre marking and numerical marking every 5cm (from 20 to 105cm).

They have 9 lateral perforations. They are silicon-lubricated on the outside and each Charriere size can be identified by a colour code on the funnel.



Product Code	Ch/Fr	Cap tether color-code	Length (cm)	Sales Unit
AS210	10	Light grey	110	20
AS410			130	
AS212			110	
AS412	12	White	130	
AS214			110	
AS414			130	
AS216	14	Dark green	110	
AS416			130	
AS218			110	
AS418	16	Orange	130	
AS220			110	
AS420			130	

MD class: IIa, CE0459. Please read the instruction for use

Helapet Limited.

Lister House, Blackburn Road, Houghton Regis, LU5 5BQ, UK
T: +44 (0)1582 501980 F: +44 (0)1582 501981 E: sales@helapet.co.uk
helapet.co.uk



Be part of our social community
by following us on:



510/0825