



Neutraclear™

Needle-free Connector

Small, compact with a reverse split septum and sealing ring which provides an effective closed system and microbiological barrier to reduce the risk of contamination. Completely transparent with a low priming volume, Neutraclear™ is suitable for paediatric doses, precise flushing and blood sampling.

Neutral fluid displacement when disconnected minimises the risk of blood reflux and occlusions occurring when used with a catheter. This helps to save on time required for any additional treatments and potential catheter restarts.

Features and Benefits

- ✓ Swabable septum surface
- ✓ Low priming volume
- ✓ Closed membrane connection reducing risk of contamination
- ✓ Fully transparent for clear view of any residual fluids
- ✓ Neutral displacement preventing occlusion and infection
- ✓ Optional end cap for additional protection



IV0200



IV0200C

Ordering Information

IV0200

Neutraclear™ Needle-free Connector
sterile • 200 pcs/pck

IV0200C

Neutraclear™ Needle-free Connector w/Cap
sterile • 200 pcs/pck

Enquire
today for
TRIAL
samples



Technical Data

	IV0200 IV0200C
Housing material	Polyoxymethylene; Polycarbonate; Silicone; Polyethylene (Permeable End Cap)
Length	2.8cm
Residual Volume	0.04ml
Properties	Neutral Displacement (-0.004ml); Disinfectable Needle-free Valve; Latex Free; DEHP Free
Maximum injection pressure	22 bars (325 psi) at a rate of 10ml/s
Back pressure	2.6 bars (37 psi)
Flow rate	>125ml/min (average)
Maximum no. of activations	400
Luer access	Swabable luer lock connection
Packaging	Individually blister packed; 200 x blisters per pack
Sterilisation	ETO
Compliance	CE marked (MDD)
Shelf Life	35 months from date of manufacture

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

