



Vented Spike with Neutraceutical™ Connector

Ideal for needle-free venting of hazardous drugs from 20mm vials, the integrated Neutraceutical™ Bi-directional connector provides a smooth swabable surface for easy disinfection. Safe and easy to use, the NeutraClear™ connector prevents the risk of drip contamination. Low priming volume makes this compact device suitable for paediatric doses.

Features and Benefits

- ✓ Secure vial attachment with finned arms
- ✓ Safe, needle-free penetration
- ✓ 0.22µm vent filter for assured sterile venting applications
- ✓ Swabable luer valve for precise aspiration of small doses
- ✓ Spike cap to reduce risk of touch contamination

*data available on request



IV0509

Ordering Information

IV0509

Vented Spike w/ Neutraceutical™ Connector
20mm • sterile • 40pcs/pck

Enquire
today for
TRIAL
samples



Technical Data

	IV0509
Housing material	Polyoxymethylene (POM); Polycarbonate (PC), Silicone (Si); Acrylonitrile butadiene styrene (ABS)
Vent Filter	0.22µm Oleo-hydrophobic
Adapter diameter	20mm
Properties	Disinfectable Surface; Latex-free
Length	7cm
Maximum injection pressure	2 bars (29 psi) at a rate of 10ml/s
Back pressure	2.6 bars (37 psi)
Flow rate	>50ml/min
Residual volume	0.14ml
Luer access	Swabable luer lock connection
Packaging	Individually blister packed; 40 x blisters per pack
Sterilisation	ETO
Compliance	CE marked (MDD)
Shelf Life	35 months from date of manufacture

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

