

HELAPET

Filter Devices • Medical & Cleanroom Disposables

To: The Editor
Date: 29 September 2003
Subject: Reduce the risk of injury with Sterile Ampoule Breaker
From: Glenn Kenworthy,
Marketing Manager, Helapet. Ltd.
Tel: 01582 501999 e.mail: glennk@helapet.co.uk

For Immediate Release

PRESS RELEASE

Reducing Risk when breaking ampoules

Helapet are pleased to add the Sterile Ampoule Breaker to their comprehensive range of products for use in hospital pharmacies, laboratories and pharmaceutical cleanrooms. The Sterile Ampoule Breaker reduces the risk of injury to staff when breaking glass ampoules. It is made from tough, rigid polypropylene, and is bright yellow to aid recognition and visibility.

Of paramount importance is safety in use. With features such as a re-closable flap for security between uses and the facility to hold up to 200 ampoule caps, makes

"This news release is issued in accordance with clause 1.2j of the British Codes of Advertising and Sales Promotion, and therefore cannot be subject to a transaction of any kind."

this a very convenient product for use throughout the hospital, however being sterile this makes it especially convenient for use in the Aseptic suite and Pharmacy isolators.

Glenn Kenworthy, Marketing Manager at Helapet notes, "This is a great little product that helps prevent injury to all staff who regularly break glass ampoules with their hands and just a tissue or wipe. Sterilising the product enables it to be used in Aseptic and cytotoxic suites as well as pharmaceutical manufacturing sites. We are convinced it will prove very useful and are pleased to offer products that can help prevent injury".

For further information and FREE samples please contact the Helapet sales team on 01582 501980 or via email at: sales@helapet.co.uk

Ends