

STOP PRESS ** STOP PRESS **** STOP PRESS**

Company: Helapet Ltd

Contact: Ruth Hewitt – 01582 501983

Date: 14th September 2009

Title: New SteriClean[®] rotational cleaning and disinfection regime now available from Helapet

Two new sterile rotational biocides are newly available to compliment the SteriClean denatured ethanol and isopropyl alcohol sprays.

Stericlean[®] Bio + I and II are good broad spectrum biocides that when used in conjunction with SteriClean[®] alcohol sprays will effectively control microbial contamination that is commonly found in cleanroom environments. Suitable for use on all hard surfaces, both inside and outside the isolators and hatches, SteriClean Bio + I and II are fully compatible with the complete range of Stericlean alcohols and wipes.

SteriClean[®] Bio + are blended with water for injection, filtered to 0.2 micron and filled in a class 10,000 cleanroom. Ready mixed and validated for maximum effect, they are supplied individually double bagged and sterilised by gamma irradiation at no less than 25 kGy.

The combined effect of the non-venting trigger and the bag in a bottle closed system, ensures the sterility of the contents throughout use and independent validation data has demonstrated that the sprays remain sterile for up to three months from opening.

Use the SteriClean[®] range of alcohols, biocides and neutral detergent to effectively control contamination in critical environments.

For more information and to try SteriClean[®] cleaning in your unit, please call our sales team on **0800 0328 428** or visit our website www.helapet.co.uk

244 words

Ends

