

STOP PRESS **** STOP PRESS **** STOP PRESS **** STOP PRESS **** STOP PRESS ****

Contact: sales@helapet.co.uk

Date: February 2017

Helapet unveil new products for 2017

Female Luer Lock Connectors / Tamper-evident IV Syringe Cap

As part of our growing portfolio of sterile syringe caps and transfer devices, Helapet Ltd have today unveiled new sterile consumable products for 2017.

[Female Luer Lock Connectors for luer lock syringes](#) facilitate safe, contaminant-free transfer of fluids from a master syringe into smaller doses. Manufactured using latex-free and non-pyrogenic materials, our double ended Luer Lock Connectors can be securely applied to male luer lock syringes preventing the risk of leaking and possible contamination.

Technicians now have more control in the sterile preparation of syringes, with a choice of Luer Lock Connectors to suit individual application. The Female Luer Lock Connector has an ergonomic, rounded design with a spacious pathway to facilitate rapid and efficient transfer of fluids between syringes. Our unique Luer Lock Micro Channel Connector includes a narrow bore pathway for precise transfer of smaller doses, leaving minimal residual waste between syringes. Delivering valuable accuracy, the Luer Lock Micro Channel Connector is perfect for paediatric medicines, small syringe sizes and high value medications.

Both Luer Lock Connectors are supplied individually blister packed making them ideal for disinfection and non-touch application to the syringe.

[Tamper-evident IV Syringe Caps for luer lock syringes](#) are a practical, cost-effective solution for the protection of syringes once they have left the sterile preparation area. An easy-to-use alternative to sticky tamper-evident tape or shrink wrap, Tamper Evident IV Syringe Caps maintain the product integrity of the syringe and allows clinical professionals to visually identify any evidence of tampering before administration.

Manufactured from latex-free materials, each cap can be applied to a luer lock syringe using the specially designed sterile tray packaging. Each sterile tray contains ten caps that are individually secured with a simple twist of the syringe, negating the risk of touch contamination and enhancing aseptic technique. Simple and easy-to-use, Helapet's new Tamper Evident IV Syringe Caps offer excellent value in maintaining drug safety.

For more information or to place an order, contact Helapet on **0800 0328 428** or email sales@helapet.co.uk. FREE samples of all our new products are now available on request.



Luer Lock Connector IV5712



Luer Lock Micro Channel Connector IV5713



Tamper Evident IV Syringe Cap IV5750R