

STOP PRESS ** STOP PRESS **** STOP PRESS **** STOP PRESS **** STOP PRESS**

Company: Helapet Ltd

Contact: Chris Steng

Date: August 2013

Title: **Your extra margin of safety with new
NeedleSafe II™ exclusively from Helapet**

In the vital drive to reduce needle-stick injury, Helapet is delighted to introduce the [NeedleSafe II™](#) as the safe, easy-to-use solution for the uncapping and re-capping of needles.

Offering superior protection, NeedleSafe II™ allows customers a simple single-handed method of uncapping and re-capping needles across a universal array of healthcare applications. Small, lightweight and highly versatile, this innovative medical device can eliminate the risk of needle-stick exposure without taking up valuable workspace.

Available exclusively to Helapet in the UK and Ireland, NeedleSafe II™ has been designed to be totally re-usable, significantly reducing medical waste and associated disposal costs. An internal gripper core allows NeedleSafe II™ to fit with several needle types including IV Catheters and Phlebotomy sets, providing essential protection whenever and wherever required.

The vast applications of NeedleSafe II™ are complemented with a detachable mounting bracket that can securely fasten the assembly horizontally or vertically onto any surface. Users can easily detach NeedleSafe II™ and clean using a number of techniques including autoclaving, ethylene oxide gas or simple alcohol disinfection.

NeedleSafe II™ has already achieved high global recognition, having been adopted to support needle-stick injury reduction in several healthcare departments. To see the benefits yourself simply call our UK freephone **0800 0328 428** or email sales@helapet.co.uk to request your free demonstration.

