

STOP PRESS ** STOP PRESS **** STOP PRESS **** STOP PRESS**

Company: Helapet Ltd

Contact: Ruth Hewitt - 01582 501983

Date: 8th September 2009

Title: **Improved specification for aseptically packaged disposables from Helapet**

Helapet have upgraded their facility by investing in a state of the art bespoke vacuum packaging machine to meet the increased demands for sterile packaging within controlled environments.

The machines are programmed per product to gain maximum efficiency and security to produce quality packaged disposables suitable for use in aseptic areas and isolator cabinets.

The machines remove the air from each pack during the packaging process ensuring that microbial particulates are extracted before sealing securely with a validated heat seal. Sterile products are then sterilised by either gamma irradiation or ETO gassing processes.

Products that now benefit from the new cleaner and more secure packaging process are, our range of PharmaPack[®] bulk packed syringes, our individually sterile wrapped sharps bins and our pharmacy prep trays.

For more information or samples of our sterile disposables for controlled environments, please call **0800 0328 428** or visit our website

www.helapet.co.uk

183 words

Ends