

PRESS RELEASE ** PRESS RELEASE **** PRESS RELEASE**

Simultaneous filtration and sterile venting in one compact device

When a needle or spike is inserted into a vial, hazardous aerosols can be released that are both harmful to operators and can cross contaminate other products in the aseptic area or isolator cabinet.

The **NEW** Rowspike Filter will protect operators and avoid cross contamination and is an ideal choice for reconstituting cytotoxics.

The Rowspike incorporates a 5 micron liquid filter and a 0.1 micron PTFE air filter to simultaneously filter a solution and sterile vent whilst simply drawing the solution into a syringe. Convenient and easy to use, your syringe attaches directly to the top of the device with a luer lock offering security whilst handling cytotoxics.

For more information on the Rowspike and other sterile venting devices, contact Helapet on **0800 0328 428** for **FREE SAMPLES**.

