

Helapet looks West for new Needleless Transfer Systems

Continuing to expand on our specialty range of syringe venting, filtration and transfer systems, Helapet are delighted to secure a partnership with West Pharmaceutical Services based in the United States.

Recognising the increasing use of needle-free systems in sterile and non-sterile clinical activities, Helapet are offering a quality range of sterile Vial Adapters and the innovative Vial2Bag® IV transfer device.

Helapet's new range of Vial adapters are designed to securely attach to any standard 20mm vial, and facilitate the optimal transfer of medical solutions for reconstitution, or transferring solutions to an IV Bag. Using a slim plastic spike for best seal penetration, the Vial Adapters can be easily applied in confined spaces including isolators and laminar flow cabinets.

The Swabable Vial Adapter is supplied with a Halkey Roberts® connector, which can maintain the sterile integrity of medicinal contents during repeated access via a needle-less luer syringe. A secure drip free system and swabable contact surface, allows the swabable vial adapter to greatly reduce the risk of cross contamination.

A new and innovative solution for transferring medical solutions from a vial or syringe into any IV bag, the Vial2Bag® is perfect for the accurate application of dosages. Using a Halkey Roberts® drip free fastening, the Vial2Bag® can significantly help reduce cross contamination with its swabable surface, whilst the inclusion of a plastic spike for IV port penetration eliminates the risk of sharps injury. The Vial2Bag® IV transfer system is ideal for use in sterile and non-sterile departments including Wards and Aseptic Production Units (APU's), supporting all IV reconstitution and administration procedures.

To learn more about Helapet's exciting new range of needleless transfer systems, please contact us on **0800 0328 428** or email sales@helapet.co.uk



Vial2Bag® IV Transfer System



Swabable Vial Adapter