

Protect products and staff from harmful aerosols

When a needle or spike is inserted into a vial, hazardous aerosols can be released that are both harmful to operators and can cross contaminate other products in the aseptic area.

Helapet offer a range of sterile venting devices to protect operators and avoid cross contamination. PharmaVent devices have a 0.2 micron filter for sterile venting and a PTFE membrane to prevent aerosols from escaping at the same time. The devices are CE marked, sterile and individually wrapped for convenient use in aseptic environments.

If a spike is preferred, they have the advantage of being suitable for multi-dose applications and offer excellent flow rates. The Taky Spike will not only sterile vent through the 0.2 micron filter, but will filter the solution through a 5 micron filter simultaneously.

For more information on sterile venting devices, contact Helapet on 0800 0328 428 for FREE SAMPLES.

