

PRESS RELEASE ** PRESS RELEASE **** PRESS RELEASE**

SteriClean Prep Pads – perfect for vials and caps

Impregnated with 70% sterilised alcohol in an individual foil sachet, SteriClean Prep Pads are exactly the right size to wipe the tops of vials or caps prior to insertion with needles or spikes.

Economic and convenient to use, the wipes are for single use and facilitate asepsis with minimal wastage.

Each 65 x 30mm wipe is impregnated with a preservative solution of 70% v/v isopropyl alcohol in water and are gamma irradiated to not less than 25KgY.

For more information or **FREE SAMPLES** of the prep pads or other impregnated wipes in the SteriClean range, please contact Helapet on **0800 0328 428**