

STOP PRESS **** STOP PRESS **** STOP PRESS **** STOP PRESS **** STOP PRESS

Company: Helapet Ltd
Contact: sales@helapet.co.uk
Date: November 2015

SteriClean® FloWrap Impregnated Wipes

SteriClean® FloWrap wipes offer exceptional cleaning and sanitising performance for cleanrooms and isolators.

Impregnated with high performance Isopropyl alcohol IPA or Denatured Ethanol (IMS) in a Class 10,000 rated cleanroom, Helapet FloWrap wipes are designed to be tough on contamination without leaving undesirable residues.

Each low linting wipe comprises of a tough, highly absorbent 68 gsm polyester/cellulose material, allowing quick and rapid disinfection of surfaces without the risk of shedding which could compromise the sterile environment.

Ready folded in 10 and 100 wipe multi-packs and individually accessed from the paper-free pouch, allowing staff to easily manage their usage.

Produced in an ISO Class 5 environment, each FloWrap pack is double bagged and gamma sterilised to no less than 25KgY supporting aseptic transfer into critically controlled environments. 'Peel and reseal' labels allow each pack to be sealed between wipes, avoiding the risk of unused wipes drying.

Users are left with only a thin piece of polythene to dispose of once the pack is fully used, ensuring less waste to dispose of and maintaining valuable work space. This makes Helapet FloWrap wipes an excellent option for confined environments such as isolators and an ideal alternative to bulky tub wipes.

Whether you require our FloWrap impregnated wipes, dry wipes or individual alcohol wipes, you can rest assured we have a solution to meet your needs.

Order your free samples of our FloWrap range by contacting us on **+44 (0)800 0328 428** or email sales@helapet.co.uk.

